



KAYDIA
WELLNESS

User Manual

Natural Tension and Headache
Relief Acupressure Hand Clip

Features:

- Kaydia provides positive pressure on the HeGu Li-4 acupressure point for natural and drug-free tension and headache relief
- Extra-large copper acupressure contact points promote blood circulation and are anti-bacterial
- Precise, consistent pressure for ease of use
- Ambidextrous & one-size-fits-most ergonomic design
- Comfortable, low-profile design for those with an active lifestyle

Care & Maintenance:

- There is a protective coating on the copper acupressure points to keep it in good condition, it will soon wear off while being used so that you can benefit from the copper properties.
- Oxidation of copper is a natural occurrence, do not be alarmed when the metal surfaces starts to change in color, it can be cleaned easily.
- Regularly wash your Kaydia Clip with soap and water.
- When oxidation/tarnishes start to show, take a soft cloth to buff the copper components with either vinegar, lemon/lime juice and or baking soda. Rinse with water afterwards.

Contents:


- 1 x Kaydia acupressure hand clip

Using Instructions:


1. Hold out your hand




2. Move your thumb to rest along your index finger. At the end of the crease between your thumb and index finger is your HeGu Li-4 point.



3. Pinch Kaydia Hand Clip onto Li4 point on hand. Make sure you can read the words "RIGHT" on the top of your right hand for correct position.



4. Apply the same for your left hand. The words "LEFT" on the top of your left hand for correct position.



WARNING

- Not suitable for expectant mothers
- Do not wear for over 30 minutes at a time

Size: 180 x 180mm
105銅版,對折4次

